



Trans, The Highest Quality  
Products For Life Sciences

## TransStart® Taq DNA Polymerase

Cat. No. AP141

Concentration 2.5 units/μl

Storage at -20°C for two years

### Description

*TransStart® Taq* DNA Polymerase is a hot start Taq DNA polymerase containing *Taq* DNA polymerase and two proprietary DNA binding proteins. At room temperature, one binding protein binds to double-strand DNA template and another binding protein binds to primer. These unique formulations effectively neutralize the DNA polymerase activity at room temperature. Blocking proteins are released from templates and primers during the initial denaturation. This double blocking method has higher efficiency than antibody based, or chemically modified hot start PCR.

### Highlights

- *TransStart® Taq* DNA Polymerase offers 18-fold fidelity as compared to *EasyTaq®* DNA Polymerase.
- Extension rate is about 1-2 kb/min.
- Template-independent “A” can be generated at the 3’ end of the PCR product. PCR products can be directly cloned into *pEASY®*-T vectors.
- Reduced nonspecific amplification and primer dimer formation.
- Different from *Taq* antibody, no risk of contamination from mammalian DNA.
- Different from chemical modification, long denaturing step is not needed.
- Amplification of genomic DNA fragment up to 15 kb.

### Applications

- Complex templates
- GC/AT-rich templates
- Multiplex PCR
- High yield PCR

### Unit Definition

One unit of *TransStart® Taq* DNA Polymerase incorporates 10 nmol of deoxyribonucleotide into acid-precipitable material in 30 minutes at 74°C.

### Quality Control

*TransStart® Taq* DNA Polymerase has passed the following quality control assays: functional absence of double- and single-strand endonuclease activity; >99% homogeneous measured by SDS-PAGE. Each batch of *TransStart® Taq* DNA Polymerase has been assayed for amplification efficiency to amplify p53 gene from 10 ng of human genomic DNA.

### Storage Buffer

20 mM Tris-HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, 100 mM KCl, 50% glycerol, stabilizers

### 10×*TransStart® Taq* Buffer with 20 mM MgSO<sub>4</sub>

500 mM Tris-HCl (pH 9.0), 200 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 20 mM MgSO<sub>4</sub>, 10% glycerol, others

### GC Enhancer

For better amplification of GC/AT-rich or complex templates, we recommend adding GC enhancer into PCR reaction. GC enhancer is provided at 10× concentration and can be used at 0.5×-5× concentration.

### Kit Contents

Component	AP141-01/11	AP141-02/12	AP141-03/13
<i>TransStart</i> <sup>®</sup> <i>Taq</i> DNA Polymerase	250 U×1	500 U×1	500 U×6
10× <i>TransStart</i> <sup>®</sup> <i>Taq</i> Buffer	1.2 ml	1.2 ml×2	1.2 ml×12
2.5 mM dNTPs	-/800 µl×1	-/800 µl×2	-/1.2 ml×8
10×GC Enhancer	200 µl×1	400 µl×1	1 ml×1
6×DNA Loading Buffer	500 µl×1	1 ml×1	1 ml×2

### Reaction Components

Component	Volume	Final Concentration
Template DNA	Variable	as required
Forward Primer (10 µM)	1 µl	0.2 µM
Reverse Primer (10 µM)	1 µl	0.2 µM
10× <i>TransStart</i> <sup>®</sup> <i>Taq</i> Buffer	5 µl	1×
2.5 mM dNTPs	4 µl	0.2 mM
<i>TransStart</i> <sup>®</sup> <i>Taq</i> DNA Polymerase	0.5-1 µl	1.25-2.5 units
ddH <sub>2</sub> O	Variable	-
Total volume	50 µl	-

### Thermal cycling conditions

94°C	2-5 min	} 30-35 cycles
94°C	30 sec	
50-60°C	30 sec	
72°C	1-2 kb/min	
72°C	5-10 min	

### Notes

- A final concentration of 2 mM MgSO<sub>4</sub> is sufficient for most targets amplification. For some targets, more Mg<sup>2+</sup> may be required.
- For optimal results, we recommend to use the 100 mM MgSO<sub>4</sub> stock to prepare a titration from 2 mM to 4 mM (final concentration) in 0.25 mM increments.
- 0.5 µl (2.5 units) enzyme is enough for per 50 µl reaction. For better amplification, up to 1 µl (5 units) enzyme can be used.
- For amplification of GC/AT-rich templates and complex templates, we suggest to use GC Enhancer.

FOR RESEARCH USE ONLY